

InBody Compatible Printer List

Manufacturer	Model	Type	Language/ Setting	Interface	S10	230	270	320	370	520	570	720	770	S10
HP	Officejet 6000	Ink	PCL Compatible	USB	○	○	○	N/C	○	○	○	○	○	○
	Officejet 6100		PCL Compatible	USB	○	○	○	N/C	○	○	○	○	○	○
	Officejet 8100		PCL Compatible	USB	○	○	○	N/C	○	○	○	○	○	○
	Mobile 470		PCL Compatible	USB	○	○	○	N/C	○	○	○	○	○	○
	Moible 100		PCL Compatible	USB	○	○	○	N/C	○	○	○	○	○	○
SAMSUNG	ML-2165W	Laser	SPL 2009/ SPL 2009-600	USB	○	○	○	N/C	○	○	○	○	○	○
	ML-2955ND		SPL 2009/ PCL Compatible	USB	○	○	○	N/C	○	○	○	○	○	○
	M2020(SL-M2020)		SPL 2009/ SPL 2009-600	USB	○	○	○	N/C	○	○	○	○	○	○
	ML-2525		SPL 2009	USB	○	○	○	N/C	○	○	○	○	○	○
	ML-2855ND		SPL 2009	USB	○	○	○	N/C	○	○	○	○	○	○
BROTHER	HL-2270DW	Laser	PCL Compatible	USB	○	○	○	○*	○	○	○	○	○	

* Requires PT-300 to convert parallel to USB

****NOTE: Please make sure your InBody has the language setting available prior to purchasing the printer*****

InBody